

Pharmacy Refrigerators *PRQ 7000 Series*



Pharmacy Refrigerators PRQ 7000 Series

Pharmacy Refrigerators PRQ 7000 series designed with good heat preservation effect, uniform internal temperature and USB port for temperature data downloading standard on upright models. It will automatically switch to energy-saving mode, under low ambient temperature, while single system runs ensuring normal operation. Under high ambient temperature, it switches to performance mode and both systems work at same time.

➤ **Design**

Designed with 7-inch touch LCD screen, all operations are operated through the screen, wide range operating voltage system from 187V to 242V allows units installed in areas with poor voltage condition. It is capable of producing two types of alarms – audible buzzer and visible flashing light.

➤ **Security System**

Equipped with safety alarm such as high/low temperature, temperature controller failure system, power failure system, abnormal door opening and closing alarm, condenser high temperature alarm, ambient temperature high and low alarm, system abnormal alarm and remote alarm interface. Also equipped with electronic lock which can unlock and open the door by fingerprint, swipe card, and account login.

➤ **Optional Accessories**

Chart Recorder

➤ **Applications**

Used in life science, healthcare, research, quality control and industrial laboratories for cold storage of drugs, vaccines, biological reagents, medicines, temperature sensitive samples, chemicals, pharmaceutical products, etc.

Specifications

Model No.	PRQ 7000	PRQ 7001
Capacity	515L	1105L
Temperature Range	2°C ~ 8°C	2°C ~ 8°C
Specialty	Completely shielded from light, specially designed for items that need to be shielded from light	
Internal Dimensions (W*D*H)	595*641*1475mm	1225*641*1562mm
External Dimensions (W*D*H)	695*794*1993mm	1325*794*1989mm
Power	350W	550W
Net Weight	155kg	200kg

Model No.	PRQ 7002	PRQ 7003
Capacity	516L	1106L
Temperature Range	2°C ~ 8°C	2°C ~ 8°C
Specialty	Visually monitor the status of internal items, and heated glass door to prevent condensation	
Internal Dimensions (W*D*H)	595*641*1475mm	1225*641*1562mm
External Dimensions (W*D*H)	695*794*1993mm	1325*794*1989mm
Power	350W	550W
Net Weight	155kg	200kg