

Under-Counter Freezers UCF 5660 Series



Under-Counter Freezers

Under-Counter Freezers UCF 5660 Series

Under-Counter Freezers UCF5660 series is equipped with LED digital display and embraco compressors. The temperature range of these freezers is -14°C to -18°C and comes in a wide range of different no. of shelves and doors. These are easy to operate, easy to maintain and can exhibit better efficiency.

Features

- ◆ R404a is the refrigerant which is free from CFC and foam
- ◆ Electronic thermostat and LED digital display for accuracy and easy reading
- ◆ 100% copper condenser and evaporator
- ◆ 20FT/40HQ are 20ft and 40ft refrigerate containers
- ◆ Auto defrosting and auto evaporation or condensation of water by means of hot pipe in stainless steel box means no floor drain required
- ◆ The internal and external structure is of stainless steel with rounded edges for easy sanitation
- ◆ Insulation 50mm thickness for better energy efficiency
- ◆ Self-closing door by means of spring-loaded hinges ensures fast recovery and energy saving
- ◆ Magnetic-Gasket with high quality neoprene rubber can be replaced without any tools Heating element concealed in door frame to eliminate any condensation
- ◆ Shelf and shelf support are designed to accommodate 1 number of GN 1/1 pan size
- ◆ The doors are structured in offset profile for easy gripping and opening as well as easy maintenance
- ◆ Ventilated refrigeration compartment gives better air circulation, anticorrosion-treated evaporator

Under-Counter Freezers

Applications _____

Used mainly in General Lab storage, blood storage, hazardous material storage, Schools & Administrative office, hyper markets and others.

Specification _____

Model No.	UCF 5660	UCF 5661	UCF 5662
Capacity	260L	260L	260L
Temperature Range	-14°C--18°C	-14°C--18°C	-14°C--18°C
No. of Doors	2	2	2
No. of 20FT/40HQ	24/48	24/48	24/48
No. of Drawers	0	2	1
Adjustable Shelves	2	2	2
Product Dimensions	1360*700*850 mm	1360*700*850 mm	1360*700*850 mm
Package Dimensions	1420*740*910 mm	1420*740*910 mm	1420*740*910 mm
Net Weight	110 kg	110 kg	111 kg
Gross Weight	135 kg	123 kg	124 kg
Power Supply	220V/50Hz or 110V/60Hz	220V/50Hz or 110V/60Hz	220V/50Hz or 110V/60Hz

Under-Counter Freezers

Model No.	UCF 5663	UCF 5664	UCF 5665
Capacity	260L	260L	260L
Temperature Range	-14°C~-18°C	-14°C~-18°C	-14°C~-18°C
No. of Doors	2	2	2
No. of 20FT/40HQ	24/48	24/48	24/48
No. of Drawers	3	4	6
Adjustable Shelves	2	2	2
Product Dimensions	1360*700*850 mm	1360*700*850 mm	1360*700*850 mm
Package Dimensions	1420*740*910 mm	1420*740*910 mm	1420*740*910 mm
Net Weight	113 kg	123 kg	129 kg
Gross Weight	126 kg	136 kg	141 kg
Power Supply	220V/50Hz or 110V/60Hz	220V/50Hz or 110V/60Hz	220V/50Hz or 110V/60Hz

Optional _____

- ◆ Set drawers instead of doors
- ◆ Castors
- ◆ Stainless steel shelf or plastic-coated shelf is available